

AOHUA



VETERINARY

Endoscopic Accessories

Aohua Endoscopy Co.,Ltd.

E-mail : vets@aohua.com

Website: www.aohua.com/vet.en



China

No.66, Lane133, Guangzhong Road,
Minhang District, Shanghai

Tel : 86-21-6768 1018

Germany

Fichtenstrasse 27,
Brunnthal, Munich

Tel :+49 (0) 81 04/8908-10

USA

20 N Front St
Bally, PA 19503

Tel : 888-238-5214

ISO 9001/ISO 13485 Certified
v2.0_2209

CONTENTS



MAKE ENDOSCOPY

DIAGNOSIS AND TREATMENT

MORE ACCESSIBLE AND CONVENIENT

NITENTS

01

Biopsy Forceps

01

02

Grasping Forceps

03

03

Accessories for
HF Surgery

06

04

Cleaning Brush

07

05

Veterinary Stent

08

06

Other Accessories

09

07

Accessories Kit

10







Biopsy Forceps









Biopsy Forceps, Sterile, Pack of 10

Disposable

	Type	Model NO.	Max insertion Portion Diameter & Working Length (mm)	Working Channel (mm)
	Oval Cup without Spike Uncoated	JRQ-Y-1816-PA	φ1.8x1600	≥2.0
		JRQ-Y-1823-PA	φ1.8x2300	≥2.0
		JRQ-Y-2323-PA	φ2.3x2300	≥2.8
		JRQ-Y-2340-PA	φ2.3x4000	≥2.8
	Oval Cup with Spike Uncoated	JRQ-Z-2323-PA	φ2.3x2300	≥2.8
		JRQ-Z-2340-PA	φ2.3x4000	≥2.8
	Alligator Cup Uncoated	JRQ-Y-1816-EA	φ1.8x1600	≥2.0
		JRQ-Y-1823-EA	φ1.8x2300	≥2.0
		JRQ-Y-2323-EA	φ2.3x2300	≥2.8
		JRQ-Y-2340-EA	φ2.3x4000	≥2.8
	Alligator Cup with Spike Uncoated	JRQ-Z-2323-EA	φ2.3x2300	≥2.8
		JRQ-Z-2340-EA	φ2.3x4000	≥2.8


Unsterile, Pack of 1

Reusable

Type	Model NO.	Max insertion Portion Diameter & Working Length(mm)	Working Channel (mm)
 Long Oval Cup not Fenestrated	AH100-10-160	φ1.0x1600	≥1.2
 Long Oval Cup Serrated	AH-118-18-160	φ1.8x1600	≥2.0
	AH-118-18-230	φ1.8x2300	≥2.0
	AH-118-23-230	φ2.3x2300	≥2.8
	AH-118-23-400	φ2.3x4000	≥2.8
	AH-118-34-400	φ3.4x4000	≥3.7
 Long Oval Cup Fenestrated	AH-100-18-160	φ1.8x1600	≥2.0
	AH-100-18-230	φ1.8x2300	≥2.0
	AH-100-23-230	φ2.3x2300	≥2.8
	AH-100-23-400	φ2.3x4000	≥2.8
	AH-100-30-400	φ3.0x4000	≥3.7
 Long Oval Cup Fenestrated with Spike	AH-101-18-160	φ1.8x1600	≥2.0
	AH-101-18-230	φ1.8x2300	≥2.0
	AH-101-23-230	φ2.3x2300	≥2.8
	AH-101-23-400	φ2.3x4000	≥2.8
	AH-101-30-400	φ3.0x4000	≥3.7
 Oval Cup with fine teeth	AH-136-18-160	φ1.8x1600	≥2.0
	AH-136-18-230	φ1.8x2300	≥2.0
	AH-136-23-230	φ2.3x2300	≥2.8
	AH-136-23-400	φ2.3x4000	≥2.8
 Round Cup fenestrated	AH-106-18-160	φ1.8x1600	≥2.0
	AH-106-18-230	φ1.8x2300	≥2.0
	AH-106-23-230	φ2.3x2300	≥2.8
	AH-106-23-400	φ2.3x4000	≥2.8
	AH-106-30-400	φ3.0x4000	≥3.7

Cytology brush, Sterile, Pack of 10

Disposable





Type	Model NO.	Max insertion Portion Diameter & Working Length(mm)	Working Channel (mm)
 Straight With sheath	JRS-A φ1.2×1200-PN	φ1.1x1200	≥1.2
	JRS-A φ1.8×1200-PN	φ1.8x1200	≥2.0
	JRS-A φ1.8×2300-PN	φ1.8x2300	≥2.0
	JRS-A φ1.8×4000-PN	φ1.8x4000	≥2.0

Grasping Forceps




Grasping Forceps, Sterile, Pack of 10

Disposable

Type	Model NO.	Max insertion Portion Diameter & Working Length(mm)	with opening width(mm)	Working Channel (mm)
 Rat teeth and Alligator type	JRY-QF-1821	φ1.8x2100	≥8.75	≥2.0
	JRY-QE-2321	φ2.3x2100	≥7	≥2.8
 Alligator type	JRY-QE-2340	φ2.3x4000	≥7	≥2.8
	JRY-QD-2321	φ2.3x2100	≥7	≥2.8
 Rat teeth type	JRY-QD-2340	φ2.3x4000	≥7	≥2.8
	JRY-QT-2321	φ2.3x2100	≥14.4	≥2.8
 Pelican Type	JRY-QT-2340	φ2.3x4000	≥14.4	≥2.8







Stone Basket, Sterile, Pack of 5

Disposable

Type	Model NO.	Max insertion Portion Diameter & Working Length(mm)	Working Channel (mm)
 Diamond shape	JRY-WD-1818-25	φ1.8x1800	≥2.0
	JRY-WD-1823-25	φ1.8x2300	≥2.0
	JRY-WD-2423-35	φ2.4x2300	≥2.8
	JRY-WD-2440-35	φ2.4x4000	≥2.8






Grasping Forceps, Unsterile, Pack of 1

Reusable

Type	Model NO.	Max insertion Portion Diameter & Working Length(mm)	Working Channel (mm)
 Alligator Forceps	AH-109-10-160	φ1.0x1600	≥1.2
	AH-109-18-160	φ1.8x1600	≥2.0
	AH-109-18-230	φ1.8x2300	≥2.0
	AH-109-23-230	φ2.3x2300	≥2.8
	AH-109-23-400	φ2.3x4000	≥2.8
	AH-109-30-400	φ3.0x4000	≥3.7
 Rat Tooth 2x1 teeth	AH-108-10-160	φ1.0x1600	≥1.2
	AH-108-18-160	φ1.8x1600	≥2.0
	AH-108-18-230	φ1.8x2300	≥2.0
	AH-108-23-230	φ2.3x2300	≥2.8
	AH-108-23-400	φ2.3x4000	≥2.8
	AH-108-30-400	φ3.0x4000	≥3.7
 3 prongs Grasper	AH-112-10-160	φ1.0x1600	≥1.2
	AH-112-18-160	φ1.8x1600	≥2.0
	AH-112-18-230	φ1.8x2300	≥2.0
	AH-112-23-230	φ2.3x2300	≥2.8
	AH-112-23-400	φ2.3x4000	≥2.8
 3 prongs Retriever round	AH-114-18-160	φ1.8x1600	≥2.0
	AH-114-18-230	φ1.8x2300	≥2.0
	AH-114-23-230	φ2.3x2300	≥2.8
	AH-114-23-400	φ2.3x4000	≥2.8
 4 prongs Retriever sharp	AH-115-18-160	φ1.8x1600	≥2.0
	AH-115-18-230	φ1.8x2300	≥2.0
	AH-115-23-230	φ2.3x2300	≥2.8
	AH-115-23-400	φ2.3x4000	≥2.8
 Foregin Body Retriever 2x2 teeth	AH-110-23-230	φ2.3x2300	≥2.8
	AH-110-23-400	φ2.3x4000	≥2.8

Grasping Forceps, Unsterile, Pack of 1

Reusable


Type	Model NO.	Max insertion Portion Diameter & Working Length (mm)	Working Channel (mm)
 Foregin Body Retriever 2x1 teeth	AH-111-18-160	φ1.8x1600	≥2.0
	AH-111-18-230	φ1.8x2300	≥2.0
	AH-111-23-230	φ2.3x2300	≥2.8
	AH-111-23-400	φ2.3x4000	≥2.8
 Grasping Forceps Serrated	AH-113-18-160	φ1.8x1600	≥2.0
	AH-113-18-230	φ1.8x2300	≥2.0
	AH-113-23-230	φ2.3x2300	≥2.8
	AH-113-23-400	φ2.3x4000	≥2.8
 Grasping Forceps for nutrition probes	AH-131-18-160	φ1.8x1600	≥2.0
	AH-131-18-230	φ1.8x2300	≥2.0
	AH-131-25-230	φ2.5x2300	≥2.8
	AH-131-25-400	φ2.5x4000	≥2.8
 Polycatch - Forceps	AH-132-25-230	φ2.5x2300	≥2.8
 Grasping Forceps Rat Tooth 2:1 teeth cross serrated	AH-148-23-230	φ2.3x2300 Opening of jaws: 14,9 mm Length of jaws: 8,4 mm	≥2.8
	AH-148-23-400	φ2.3x4000 Opening of jaws: 14,9 mm Length of jaws: 8,4 mm	≥2.8

Accessories for HF Surgery




Snare, Sterile, Pack of 5

Disposable

Type	Model NO.	Max insertion Portion Diameter & Working Length(mm)	loop width(mm)	Working Channel (mm)
 Oval	JRD-QA-1823-25	$\phi 1.8 \times 2300$	25	≥ 2.0
	JRD-QA-2423-35	$\phi 2.4 \times 2300$	35	≥ 2.8
	JRD-QA-2440-35	$\phi 2.4 \times 4000$	35	≥ 2.8


Injection Needle, Sterile, Pack of 5

Disposable


Type	Model NO.	Max insertion Portion Diameter & Working Length(mm)	Working Channel (mm)	Needle's diameter	Needle's length (mm)
 Injection Needle	JRZ-1816-2504	$\phi 1.8 \times 1600$	≥ 2.0	25G	4
	JRZ-1823-2504	$\phi 1.8 \times 2300$	≥ 2.0	25G	4
	JRZ-2323-2305	$\phi 2.3 \times 2300$	≥ 2.8	23G	5
	JRZ-2340-2305	$\phi 2.3 \times 4000$	≥ 2.8	23G	5

Aspiration catheters, Sterile, Pack of 10

Disposable


Type	Model NO.	Max insertion Portion Diameter & Working Length(mm)	Working Channel (mm)
 Aspiration Catheters	ACT120S	$\phi 1.0 \times 1200$	≥ 1.2
	ACT200S	$\phi 1.8 \times 2000$	≥ 2.0
	ACT230S	$\phi 2.3 \times 2300$	≥ 2.8
	ACT400S	$\phi 2.3 \times 4000$	≥ 2.8

Cleaning Brush





Cleaning Brushes, Sterile, Pack of 5

Disposable

Type	Model NO.	Max insertion Portion Diameter & Working Length (mm)	Bristle Diameter (mm)	Working Channel (mm)
 <p>Double end Cleaning Brushes</p>	JRX1012-0210	φ1.0x1200	2x10	≥1.2
	JRX1812-0320	φ1.8x1200	3x20	≥2.0
	JRX1823-0320	φ1.8x2300	3x20	≥2.0
	JRX1823-0520	φ1.8x2300	5x20	≥2.8
	JRX1840-0520	φ1.8x4000	5x20	≥2.8
	JRX1840-0620	φ1.8x4000	6x20	≥3.7



Cleaning Brushes, Sterile, Pack of 1

Reusable

Type	Model NO.	Max insertion Portion Diameter & Working Length (mm)	Bristle Diameter (mm) x Length (mm)	Working Channel (mm)
 <p>Double end Cleaning Brushes</p>	JRF1823-0512	φ1.8x2300	5x12	≥2.8
	JRF1840-0512	φ1.8x4000	5x12	≥2.8
 <p>Valves Brushes</p>	JRF95-1040	φ1.0x95	10x40	/

OS

Veterinary Stent



Veterinary Tracheal Stent (uncoated straight shape), Niti-alloy Wire With Introducer

Disposable



Model NO.	Diameter(mm)	Working Length(mm)
GS-8-11-08030	8	30
GS-8-11-08050	8	50
GS-8-11-08070	8	70
GS-8-11-08080	8	80
GS-8-11-08090	8	90
GS-8-11-10030	10	30
GS-8-11-10050	10	50
GS-8-11-10060	10	60
GS-8-11-10070	10	70
GS-8-11-10080	10	80
GS-8-11-10090	10	90
GS-8-11-12040	12	40
GS-8-11-12050	12	50
GS-8-11-12060	12	60
GS-8-11-12070	12	70
GS-8-11-12090	12	90
GS-8-11-14030	14	30
GS-8-11-14040	14	40
GS-8-11-14060	14	60
GS-8-11-14070	14	70
GS-8-11-14080	14	80
GS-8-11-14100	14	100

Other Accessories

Mucus extractors, individually packed and Sterile, Pack of 5

Disposable



Type

BAL / Tracheal wash catheters can become blocked with thick mucus, collection traps ensure you retrieve the best sample to aid in your diagnosis. Single chamber mucus extractors are;

Gatroscopy Hood

Reusable




Type

In order to reduce the risk of damage to the Cardiac sphincter, Esophagus and Pharynx during foreign body retrieval our protective hood can be fitted to the majority of scopes to ensure safe extraction of a foreign body via gatroscopy. Hood Measurements: 40mm at Distal end, 8.3mm at proximal end, length of 75mm.

Guttural Pouch Probe

Reusable



Diameter(mm)	Working Length(mm)
1.8	1800
2.4	1800
2.4	2250

Accessories Kit



Disposable

Reusable

VET-3512



QTY	Type	Features	Model No.	Diameter (mm)	Working Length(mm)
<input type="checkbox"/>	Cytology Brush	Straight With sheath	JRS-A Φ 1.2 \times 1200-PN	1.2	1200
<input type="checkbox"/>	Aspiration Catheters		ACT120S	1.0	1200
<input type="checkbox"/>	Biopsy Forceps	Long Oval Cup	AH100-10-160	1.0	1600
<input type="checkbox"/>	Grasping Forceps	Alligator	AH-109-10-160	1.0	1600
<input type="checkbox"/>	Grasping Forceps	Rat Tooth 2x1 teeth	AH-108-10-160	1.0	1600
<input type="checkbox"/>	Grasping Forceps	3 prongs Grasper	AH-112-10-160	1.0	1600

VET-6010




QTY	Type	Features	Model No.	Diameter (mm)	Working Length(mm)
<input type="checkbox"/>	Biopsy Forceps	Oval Cup	JRQ-Y-1816-PA	1.8	1600
<input type="checkbox"/>	Biopsy Forceps	Alligator Cup	JRQ-Y-1816-EA	1.8	1600
<input type="checkbox"/>	Cytology Brush		JRS-A Φ 1.8 \times 1200-PN	1.8	1200
<input type="checkbox"/>	Aspiration Catheters		ACT200S	1.8	2000
<input type="checkbox"/>	Grasping Forceps	Rat teeth & Alligator type	JRY-QF-1821	1.8	2100
<input type="checkbox"/>	Stone basket	Diamond shape	JRY-WD-1818-25	1.8	1800
<input type="checkbox"/>	Snare	Oval	JRD-QA-1823-25	1.8	2300
<input type="checkbox"/>	Injection Needle		JRZ-1816-2504	1.8	1600

Accessories Kit


Disposable

Reusable

VET-8015

	QTY	Type	Features	Model No.	Diameter (mm)	Working Length(mm)
	<input type="checkbox"/>	Biopsy Forceps	Oval Cup	JRQ-Y-1823-PA	1.8	2300
	<input type="checkbox"/>	Biopsy Forceps	Alligator Cup	JRQ-Y-1823-EA	1.8	2300
	<input type="checkbox"/>	Cytology Brush		JRS-A Φ 1.8 \times 2300-PN	1.8	2300
	<input type="checkbox"/>	Grasping Forceps	Rat teeth & Alligator type	JRY-QF-1821	1.8	2100
	<input type="checkbox"/>	Stone basket	Diamond shape	JRY-WD-1823-25	1.8	2300
	<input type="checkbox"/>	Snare	Oval	JRD-QA-1823-25	1.8	2300
	<input type="checkbox"/>	Injection Needle		JRZ-1823-2504	1.8	2300

VET-8028


	QTY	Type	Features	Model No.	Diameter (mm)	Working Length(mm)
	<input type="checkbox"/>	Biopsy Forceps	Oval Cup	JRQ-Y-2323-PA	2.3	2300
	<input type="checkbox"/>	Biopsy Forceps	Oval Cup with Spike	JRQ-Z-2323-PA	2.3	2300
	<input type="checkbox"/>	Biopsy Forceps	Alligator Cup	JRQ-Y-2323-EA	2.3	2300
	<input type="checkbox"/>	Biopsy Forceps	Alligator Cup with Spike	JRQ-Z-2323-EA	2.3	2300
	<input type="checkbox"/>	Cytology Brush		JRS-A Φ 1.8 \times 2300-PN	1.8	2300
	<input type="checkbox"/>	Grasping Forceps	Rat teeth	JRY-QD-2321	2.3	2100
	<input type="checkbox"/>	Grasping Forceps	Alligator type	JRY-QE-2321	2.3	2100
	<input type="checkbox"/>	Grasping Forceps	Pelican type	JRY-QT-2321	2.3	2100
	<input type="checkbox"/>	Stone basket	Diamond shape	JRY-WD-2423-35	2.4	2300
	<input type="checkbox"/>	Snare	Oval	JRD-QA-2423-35	2.4	2300
	<input type="checkbox"/>	Injection Needle		JRZ-2323-2305	2.3	2300

Accessories Kit

Disposable

Reusable

VET-9215



QTY	Type	Features	Model No.	Diameter (mm)	Working Length(mm)
<input type="checkbox"/>	Biopsy Forceps	Oval Cup	JRQ-Y-2323-PA	2.3	2300
<input type="checkbox"/>	Biopsy Forceps	Oval Cup with Spike	JRQ-Z-2323-PA	2.3	2300
<input type="checkbox"/>	Biopsy Forceps	Alligator Cup	JRQ-Y-2323-EA	2.3	2300
<input type="checkbox"/>	Biopsy Forceps	Alligator Cup with Spike	JRQ-Z-2323-EA	2.3	2300
<input type="checkbox"/>	Cytology Brush		JRS-A Φ 1.8 \times 2300-PN	1.8	2300
<input type="checkbox"/>	Grasping Forceps	Rat teeth	JRY-QD-2321	2.3	2100
<input type="checkbox"/>	Grasping Forceps	Alligator type	JRY-QE-2321	2.3	2100
<input type="checkbox"/>	Grasping Forceps	Pelican type	JRY-QT-2321	2.3	2100
<input type="checkbox"/>	Stone basket	Diamond shape	JRY-WD-2423-35	2.4	2300
<input type="checkbox"/>	Snare	Oval	JRD-QA-2423-35	2.4	2300
<input type="checkbox"/>	Injection Needle		JRZ-2323-2305	2.3	2300
<input type="checkbox"/>	Aspiration Catheters		ACT230S	2.3	2300

VET-9833/VET-1335



QTY	Type	Features	Model No.	Diameter (mm)	Working Length(mm)
<input type="checkbox"/>	Biopsy Forceps	Oval Cup	JRQ-Y-2340-PA	2.3	4000
<input type="checkbox"/>	Biopsy Forceps	Oval Cup with Spike	JRQ-Z-2340-PA	2.3	4000
<input type="checkbox"/>	Biopsy Forceps	Alligator Cup	JRQ-Y-2340-EA	2.3	4000
<input type="checkbox"/>	Biopsy Forceps	Alligator Cup with Spike	JRQ-Z-2340-EA	2.3	4000
<input type="checkbox"/>	Cytology Brush		JRS-A Φ 1.8 \times 4000-PN	1.8	4000
<input type="checkbox"/>	Grasping Forceps	Rat teeth	JRY-QD-2340	2.3	4000
<input type="checkbox"/>	Grasping Forceps	Alligator type	JRY-QE-2340	2.3	4000
<input type="checkbox"/>	Grasping Forceps	Pelican type	JRY-QT-2340	2.3	4000
<input type="checkbox"/>	Stone basket	Diamond shape	JRY-WD-2440-35	2.4	4000
<input type="checkbox"/>	Snare	Oval	JRD-QA-2440-35	2.4	4000
<input type="checkbox"/>	Injection Needle		JRZ-2340-2305	2.3	4000
<input type="checkbox"/>	Aspiration Catheters		ACT400S	2.3	4000